

LDT 5412

Thermoelectric Temperature Controller

Product Features

High stability temperature control within 0.01°C


Hybrid P-I control loop for fast settling time

Three display modes provide easy operation

Application flexibility through user set control and limits

User selectable thermistor current

Designed to be economical and flexible, the LDT-5412 Thermoelectric Temperature Controller offers a stable, low noise bipolar output and is optimized for controlling the temperature of laser diodes and photodetectors. The instrument controls and displays thermistor resistance while delivering bipolar current to a thermoelectric module. The unit's hybrid proportional-integral control loop offers fast settling times with a typical temperature stability of 0.01°C.



Low Cost, High Stability,
Thermoelectric Temperature Control

LDT 5412

Thermoelectric Temperature Controller

Easy Operation

The intuitive front panel features a highly visible LED display, which has three display modes for easy operation: (1) Set Resistance, which displays the setpoint resistance (temperature) level, (2) Actual Resistance, which displays the actual resistance of the thermistor sensor, and (3) TE Current, which displays the drive current to the TE module.

Two visual fault indicators located on the front panel provide immediate information of fault conditions. The OPEN THERM indicator lights when the connection to the thermistor sensor is open. When this occurs, the unit's output automatically shuts off. The LIMIT indicator lights whenever the output is actively limited to the user-set current limit value.

Application Flexibility

The incorporation of user-selectable thermistor source currents of 10 μ A and 100 μ A ensures versatility over a wide range of temperatures and applications. This allows the LDT-5412 to operate over a thermistor control range of 10 Ω to 200 k Ω .

For a typical 10k Ω thermistor, this corresponds to a temperature range of -20°C to 50°C. Other temperature ranges are possible with different thermistors.

The LDT-5412 allows the user to bypass the internal set resistance function of the front panel knob and externally connect a resistor of known value. This is convenient for reproducing the same temperature quickly and accurately.

In addition, the LDT-5412 optimizes slew rate and settling time. Adjustment is easy with the rear panel GAIN control. For automated testing or to remotely compute actual temperature, the LDT-5412 also offers an analog voltage output that corresponds to the thermistor resistance.

The unit's output can be current limited anywhere within its 0-2A range by using the recessed front panel LIMIT control knob. During operation, the thermoelectric module current is unconditionally held to the set limit value, ensuring the TE module cannot be over driven.

Specifications

OUTPUT

Output Type:	Bipolar current source
Current Range:	-2 to 2A, floating
Compliance Voltage:	>2V, DC

DISPLAY

Type:	3.5 digit green LED		
	10 μ A	100 μ A	TE Current
Maximum Reading:	199.9 k Ω	19.99 k Ω	1.99A
Resolution:	0.1 k Ω	0.01 k Ω	0.01A
Accuracy (@ 25°C):	± 0.5 k Ω	± 0.05 k Ω	± 0.05 A

ACTUAL R MONITOR

	10 μ A	100 μ A
Output:	100 k Ω /V	10 k Ω /V
Accuracy:	$\pm 5\%$	$\pm 5\%$


GENERAL

AC Power:	100-120, 220-240 VAC; 50/60 Hz
Size (HxWxD):	267mm x 140mm x 66mm; 10.5" x 5.5" x 2.6"
Weight:	1.8 kg (4.0 lbs.)
Operating Temperature:	0°C to 50°C
Storage Temperature:	-40°C to 70°C
Warm-Up:	1 hour, to rated accuracy
Output Connector:	15-pin D-sub
Reference Input:	BNC connector

ORDERING INFORMATION

LDX-5412	Thermoelectric Temperature Controller (4W, includes on TS-510 Thermistor)
LDM-4405	Low Cost TO-Can Laser Diode Mount
CC-501S	TE Controller / Unterminated Interconnect Cable
CC-505S	TE Controller / Laser Diode Mount Interconnect Cable
TS-510	10k Ω Calibrated Thermistor ($\pm 0.2^\circ$ C)
TS-520	10k Ω Uncalibrated Thermistor ($\pm 1.5^\circ$ C; -20°C to 50°C)
TS-521	5k Ω Uncalibrated Thermistor ($\pm 1.5^\circ$ C; -45°C to 30°C)
TS-523	20k Ω Uncalibrated Thermistor ($\pm 1.5^\circ$ C; -10°C to 70°C)
TS-525	100k Ω Uncalibrated Thermistor ($\pm 1.5^\circ$ C; 10°C to 110°C)
UCA-350	Unipolar Heater Control Adapter

In keeping with our commitment to continuing improvement, ILX Lightwave reserves the right to change specifications without notice or liability for such changes.



With a typical drift within 0.01°C over a four-hour period, the LDT-5412 out performs similar thermoelectric temperature controllers.

ILX Lightwave

Laser Diode Instrumentation & Test Systems

P.O. Box 6310, Bozeman, MT 59771 • FAX: 406-586-9405

www.ilxlightwave.com

For information call

1-800-459-9459

International Inquiries: 406-556-2481
email: sales@ilxlightwave.com



REV04.042909