



CC306H WIRING DETAIL



- SHIELD = SHELL GND
- PR#1 < BROWN = A3 LD ANODE
RED = A3 LD ANODE
- PR#2 < WHT/BRN = A2 LD CATHODE
WHT/RED = A2 LD CATHODE
- PR#3 < ORANGE = 9 ANODE SENSE
WHT/ORG = 4 CATHODE SENSE
- PR#4 < BLACK = 8 PD ANODE
WHT/BLK = 7 PD CATHODE

CONNECTOR TYPE: 13W3 MALE

00	6/15/05	CREATED & RELEASED	SAR
REV	DATE	DESCRIPTION	CHK
REVISIONS			

UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES		FINISH			
TOLERANCES:		MATERIAL		N/A	
FRACTIONAL	± .02	ILX PART		DESCRIPTION	
DECIMAL .XX	± .005	THIS DRAWING IS THE PROPERTY OF ILX LIGHTWAVE CORPORATION. IT IS FURNISHED TO YOU FOR CONFIDENTIAL INFORMATION PURPOSES ONLY AND IS NOT TO BE DISCLOSED TO ANYONE, REPRODUCED, OR USED FOR MANUFACTURING PURPOSES WITHOUT THE WRITTEN PERMISSION OF ILX LIGHTWAVE CORPORATION.			
ANGULAR		DRAWN BY	DATE	 Telephone: (406) 586-1244 31950 E. Frontage Rd. Bozeman, Montana 59715	
HOLE TOL.		<i>S. Remington</i>	6/15/05		
DO NOT SCALE DRAWING		UPDATED BY	DATE	APPRV BY	DATE
AUTOCAD # CC306H				SCALE	TITLE
				CC306H WIRING DETAIL	
				DWG. NO. CC306H	REV 00 SHEET 1/1